
MED DIMENSIONS

Veterinary Education Catalog



Med Dimensions, LLC. || 318 Timothy Lane, Suite 2 || Ontario, NY 14519

January 2024 Edition



Overview

Founded in 2019, Med Dimensions believes that educating the next generation of veterinarians is crucial to the overall health of our pets and animals. "See one, do one, teach one" provides some real experience for students and new clinicians, but can be a slow process and limit exposure.

Med Dimensions takes real anatomy and pathology and transforms it into three dimensional solutions for any type of procedure and any specialty.

From bronchoscopy trainers to bone replication, and anything in between, we have a solution for you.

Goals

When students use a Med Dimensions model they will:

1. Reduce uncertainties in medical procedures
2. Have an increased confidence in their abilities to perform techniques and procedures
3. Experience realistic and immersive training that closely mimics real-life situations

Table of Contents

<u>Product</u>	<u>SKU</u>	<u>Page</u>
Vet Bones Canine Forelimb	MULTIPLE	5
Vet Bones Canine Hind Limb	MULTIPLE	6
Vet Bone Canine Skull Clear	MD-E-C0011Clear	7
Canine Stifle Visual Aid	MD-E-C0005E & MD-E-C0005EColor	8
Suture Pad Economy	MD-E-C1000E & MD-E-1001E	9
Suture Pad Standard	MD-E-C1000S & MD-E-1001S	10
Canine Ear	MD-E-C0011S	11
Feline Intubation Economy	MD-E-F1021E	12
Feline Intubation Standard & with NG	MD-E-F1021S & MD-E-F1031S	13
Canine Intubation Standard	MD-E-C1021S	14
Canine Brachycephalic Intubation Standard	MD-E-C1031S	15
Cockatoo Intubation Standard	MD-E-A1021S	16
Macaw Intubation Standard	MD-E-A1031S	17
Rabbit Intubation Standard	MD-E-R1021S	18
Canine Cephalic Venipuncture	MD-E-C1044S	19
Canine Jugular Access	MD-E-C1051S	20
Canine Cricothyrotomy & Tracheostomy	MD-E-C2031S	21
Canine Bronchoscopy Trainer	MD-E-C2002SM & MD-E-C2002SL	22
Canine Pericardiocentesis Economy	MD-E-C2012E	23
Canine Pericardiocentesis Standard	MD-E-C2012S	24

Table of Contents

<u>Product</u>	<u>SKU</u>	<u>Page</u>
Canine Thoracocentesis	MD-E-C2023	23
Canine Joint Injection Pelvic Limb	MD-E-C3005	24
Canine Arthroscopy Stifle	MD-E-C3015	25
Canine Vet Tech Kit	MD-E-C7001S	26
Canine Emergency Neck Access Kit	MD-E-C7001S	27
Canine Airway Management Kit	MD-E-C7011S	28
Exotic Intubation Kit	MD-E-E7021S	29

Vet Bones Canine ForeLimb



General Specifications

- Derived from CT scans of a Medium sized canine
- Cuts and drills like real bone, and can hold anchors, buttons, plates and screws

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Humerus Left	MD-E-C0014L	\$27.30
Canine Humerus Right	MD-E-C0014R	\$27.30
Canine Radius Left	MD-E-C0024L	\$27.30
Canine Radius Right	MD-E-C0024R	\$27.30
Canine Ulna Left	MD-E-C0034L	\$25.20
Canine Ulna Right	MD-E-C0034R	\$25.20
Canine Radius/Ulna Left	MD-E-C0044L	\$27.30
Canine Radius/Ulna Right	MD-E-C0044R	\$27.30
Canine Scapula Left	MD-E-C0054L	\$27.30
Canine Scapula Right	MD-E-C0054R	\$27.30

Vet Bones Canine Hind Limb



General Specifications

- Derived from CT scans of a Medium sized canine
- Cuts and drills like real bone, and can hold anchors, buttons, plates and screws

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Femur Left	MD-E-C0015L	\$27.30
Canine Femur Right	MD-E-C0015R	\$27.30
Canine Tibia Left	MD-E-C0025L	\$27.30
Canine Tibia Right	MD-E-C0025R	\$27.30
Canine Fibula Left	MD-E-C0035L	\$23.10
Canine Fibula Right	MD-E-C0035R	\$23.10
Canine Tibia/Fibula Left	MD-E-C0045L	\$27.30
Canine Tibia/Fibula Right	MD-E-C0045R	\$27.30

Vet Bones Canine Skull Clear




General Specifications

- Derived from CT scans of a Medium sized canine
- Includes Rigid Skull and C1
- Can use with dry erase markers to indicate key anatomic landmarks and enhance active education
- 120x120x230mm

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Skull Clear	MD-E-C0011Clear	\$225.75

Canine Stifle Visual Aid



General Specifications

- Derived from CT scans of a Small/medium sized canine (15-20kg and 8.5in long)
- Includes Cranial Cruciate Ligament (CCL), Caudal Cruciate Ligament (CaCL), meniscus, long digital extensor tendon, lateral collateral ligament (LCL) and medial collateral ligament (MCL)
- Ligaments will hold suture and are firmly secured within the bone
- Includes partial tibia and femur that can hold anchors, buttons, plates and screws
- 180x120x50mm Stifle based on the left hind limb

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Stifle Visual Aid	MD-E-C0005E	\$133.35
Canine Stifle Visual Aid with Color	MD-E-C0005EColor	\$147.00

Suture Pad Economy



General Specifications

- Two different materials to emulate the feeling of Dermal and Subcutaneous layers of flesh
- Embedded mesh to create appropriate “pop” feeling and improves durability
- Encased in a hard plastic case to increase longevity and stability while suturing
- Ideal for use with Monofilament Sutures
- 3 Slits 120x80mm, 6.5mm thick with 3 slits that are 100mm long
- 5 Slits 100x100mm, 6.5mm thick with 3 slits that are 40mm long, 2 that are 80mm long

<u>Component</u>	<u>SKU</u>	<u>Price</u>
3 Slit Suture Pad (Bundle of 5)	MD-E-C1000E	\$96.60
5 Slit Suture Pad (Bundle of 5)	MD-E-C1001E	\$107.10

Suture Pad Standard



General Specifications

- Three different materials to emulate the feeling of Epidermis, Dermal and Subcutaneous layers of flesh
- Embedded mesh to create appropriate “pop” feeling and improves durability
- Encased in a hard plastic case to increase longevity and stability while suturing
- Ideal for use with Monofilament Sutures
- 3 Slits 120x80mm, 6.5mm thick with 3 slits that are 100mm long
- 5 Slits 100x100mm, 6.5mm thick with 3 slits that are 40mm long, 2 that are 80mm long

<u>Component</u>	<u>SKU</u>	<u>Price</u>
3 Slit Suture Pad (Bundle of 5)	MD-E-C1000S	\$133.35
5 Slit Suture Pad (Bundle of 5)	MD-E-C1001S	\$144.90

Canine Ear




General Specifications

- Near unlimited use
- Ideal for otoscopic ear exams
- Realistic features including veins and tympanic membrane
- Comes with suction cup for easy mounting
- 150x190x90mm Full model
- 80x80x70mm Replacement ear puck

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Ear	MD-E-C0011S	\$472.50
Replacement Ear	MD-E-C0011S-1	\$393.75

Feline Intubation Economy



General Specifications

- Designed from real CT scans, and is an exact representation of real anatomy
- Near unlimited use
- Easy to manipulate with fingers, a tube, or laryngoscope
- Encased in a hard plastic case to increase longevity and stability while practicing
- Ideal for use with Size 3 ET Tube
- 95x95x130mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Feline Intubation Economy	MD-E-F1021E	\$157.50

Feline Intubation Standard

Feline Intubation Standard with NG



General Specifications

- Designed from real CT scans, with nasal cavity, teeth epiglottis, and pharynx
- NG model includes serviceable nose patch to hold suture for NG tube placement
- Replaceable jaw and tongue, and teeth on both models, NG model allows for replaceable nose
- Easy to manipulate with fingers, a tube, or laryngoscope
- Threaded to fit any standard tripod, or can be free standing
- Comes with suction cup for easy mounting
- Ideal for use with 3 ET Tube
- 140x150x90mm Full intubation model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Feline Intubation Standard	MD-E-F1021S	\$630.00
Replacement Jaw/Tongue	MD-E-F1021S-1	\$126.00
Replacement Teeth	MD-E-F1021S-2	\$26.25
Feline Intubation Standard w/ NG	MD-E-F1031S	\$735.00
Replacement Nose Puck	MD-E-F1031S-1	\$68.25

Feline Intubation Standard

Bisected



General Specifications

- Designed from real CT scans, with nasal cavity, teeth epiglottis, and pharynx
- Easy to manipulate with fingers, a tube, or laryngoscope
- Carbon copy of Feline Intubation Standard, bisected to allow for visualization of the internal structures of the feline airway
- Magnets to secure the two halves of the head
- Threaded to fit any standard tripod, or can be free standing
- Comes with suction cup for easy mounting
- Ideal for use with 3 ET Tube and 5 French NG Tube
- 140x150x90mm Full intubation model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Feline Intubation Standard Bisected	MD-E-F1041S	\$630.00

Please note that any Med Dimensions Standard Intubation Model can be bisected.
Reach out to info@med-dimensions to learn more.

Canine Intubation Standard



General Specifications

- 105x135x200mm Full model
- Designed from real CT scans, with a nasal cavity, teeth, epiglottis, and pharynx
- Replaceable jaw and tongue
- Easy to manipulate with fingers, a tube, or laryngoscope
- Threaded to fit any standard tripod, or can be free standing
- Ideal for use with 4.5 ET Tube
- Comes with suction cup for easy mounting

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Intubation Standard	MD-E-C1021S	\$945.00
Replacement Jaw/Tongue	MD-E-C1021S-1	\$189.00

Canine Brachycephalic Intubation Standard



General Specifications

- Designed from real CT scans, with a nasal cavity, epiglottis, and pharynx
- Realistic and replaceable soft palate occlusion of trachea to allow access to airway
- Replaceable jaw and tongue
- Easy to manipulate with fingers, a tube, or laryngoscope
- Threaded to fit any standard tripod, or can be free standing
- Comes with suction cup for easy mounting
- Ideal for use with 3.5 ET Tube
- 145x175x165mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Brachycephalic Intubation	MD-E-C1031S	\$997.50
Replacement Jaw/Tongue	MD-E-C1031S-1	\$194.25
Replacement Soft Palate	MD-E-C1031S-2	\$42.00

Cockatoo Intubation Standard




General Specifications

- Designed from real CT scans, with soft tongue, flexible jaw and syrinx
- Ideal for practicing intubation to the syrinx
- Ideal for use with 3.5 ET Tube
- 260x130x150mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Cockatoo Intubation Standard	MD-E-A1021S	\$367.50

Macaw Intubation Standard



General Specifications

- Designed from real CT scans, with soft tongue, flexible jaw and syrinx
- Ideal for practicing intubation to the syrinx
- Ideal for use with 5 ET Tube
- 265x155x150mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Macaw Intubation Standard	MD-E-A1031S	\$367.50

Rabbit Intubation Standard



General Specifications

- Designed from real CT scans, epiglottis, glottis, tongue and upper lip and jaw
- Ideal for practicing intubation on the difficult rabbit anatomy
- Ideal for use 3 ET Tube
- 160x130x75mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Rabbit Intubation Standard	MD-E-R1021S	\$393.75

Canine Cephalic Venipuncture



General Specifications

- Lifelike features include palpable vein, and proper positioning
- Accommodates approximately a 24 Gauge Needle
- Two types of veins allow for a variety of skill levels to practice
- High quality materials allow for repeated usage and extensive training sessions
- Replaceable skin and veins
- Comes with suction cup for easy mounting
- 70x100x300mm Full model
- 5x10x210mm Replacement Vein

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Cephalic Venipuncture	MD-E-C1044S	\$1,155.00
Replacement Skin	MD-E-C1044S-1	\$105.00
Replacement Vein Beginner 5 Count	MD-E-C1044S-2	\$131.25
Replacement Vein Standard 5 Count	MD-E-C1044S-3	\$131.25

Canine Jugular Access



General Specifications

- Realistic features of a canine's neck, including the jugular vein, skin and underlying structures
- High quality materials allow for repeated usage and extensive training sessions
- Replaceable skin and veins
- Replaceable central line cradle that holds model in dorsal recumbency
- Conversion Kit includes Jugular Access Muscle Puck, 2 Veins, and 1 Skin
- 170x370x345mm Full model
- 18x18x338mm Replacement Vein

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Jugular Access	MD-E-C1051S	\$1,260.00
Central Line Cradle	MD-E-C1051S-1	\$73.50
Replacement Skin	MD-E-C1051S-2	\$126.00
Replacement Vessels- 3 Pair	MD-E-C1051S-3	\$147.00
Conversion Kit for Cric/Trache	MD-E-C1051S-4	\$840.00

Canine Cricothyrotomy & Tracheostomy



General Specifications

- Realistic features of a canine's neck, including the Cricothyroid and Trachea with Membrane Tape
- High quality materials allow for repeated usage and extensive training sessions
- Replaceable skin, and Membrane Tape
- Conversion Kit includes Cric/Trache Muscle Puck, Airway, 6 pieces of Membrane Tape, and 1 Skin
- 170x370x345mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Cricothyrotomy & Tracheostomy Trainer	MD-E-C2031S	\$1,260.00
Replacement Membrane Tape (12 Count)	MD-E-C2031S-1	\$18.90
Replacement Skin	MD-E-C1051S-2	\$126.00
Conversion Kit for Jug Access	MD-E-C2031S-3	\$840.00

Canine Bronchoscopy Trainer



General Specifications

- Near unlimited use, and requires no assembly
- Replicates intricate anatomy structures of the respiratory system
- Flexible and responsive airway walls allow for users to navigate the bronchial tree
- Can incorporate a wide variety of simulated pathologies and anatomical anomalies
- 100x110x470mm Medium bronchoscopy trainer
- 125x160X600mm Large bronchoscopy trainer

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Bronchoscopy Medium	MD-E-C2002SM	\$735.00
Canine Bronchoscopy Large	MD-E-C2002SL	\$840.00

Canine Pericardiocentesis Economy



General Specifications

- Ultrasoundable model that can be filled with liquid to represent the pericardial sac
- Self healing material that allows for near unlimited use
- Accommodates a wide range of needle gauges
- 85x125x150mm

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Pericardiocentesis Econ	MD-E-C2012E	\$630.00

Canine Pericardiocentesis Standard



General Specifications

- Ultrasoundable model that can be filled with liquid to represent the pericardial sac
- Chest wall includes ribs in opaque material for increased realism
- Heart can be filled with water or air and can be inflated using a syringe
- Self healing material that allows for near unlimited use
- Accommodates a wide range of needle gauges
- 85x125x150mm

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Pericardiocentesis Std	MD-E-C2012S	\$787.50

Canine Thoracocentesis



General Specifications

- Replicates the thoracic cavity, from the ribs and intercostal spaces to the pleural space and lung structure
- Simulates pleural effusion, allowing identification and drainage of air from the pleural space
- High quality materials allow for repeated use and extensive practice sessions
- Embedded mesh to create appropriate “pop” feeling in the skin and improves durability
- 220x260x400mm Full model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Thoracocentesis	MD-E-C2023	\$2,835.00

Canine Joint Injection Pelvic Limb



General Specifications

- Replicates anatomy of canine pelvic limb with realistic range of motion
- Ideal for practicing the injection of stem cells in the hip and stifle joints
- Joint pucks includes realistic fluid that would be present in a live canine
- High quality materials allow for repeated usage and extensive training sessions
- Replaceable skin, and hip and stifle joints
- 145x300x420mm Full model
- 50x55x110mm Stifle Joint
- 50x57x45mm Hip Joint

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Joint Injection Pelvic Limb	MD-E-C3005	\$2100.00
Replacement Skin	MD-E-C3005-1	\$157.50
Replacement Pucks	MD-E-C3005-2	\$110.25

Canine Arthroscopy Stifle



General Specifications

- Replicates anatomy of canine stifle joint, with femur, tibia, patella, ligaments, and surrounding structures with realistic range of motion
- Stifle includes the CCL, CaCL, meniscus, long digital extensor tendon, LCL, MCL, and fat pad
- Can be used with any standard arthroscopy equipment, and has ability to run fluid through the joint
- High quality materials allow for repeated usage and extensive training sessions
- Replaceable skin and stifle joint
- 145x300x420mm Full model
- 50x55x110mm Stifle

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Arthroscopy Stifle	MD-E-C3015	\$1,260.00
Replacement Skin	MD-E-C3015-1	\$157.50
Replacement Puck	MD-E-C3015-2	\$110.25

Canine Vet Tech Kit



General Specifications

- Kit includes Canine Brachycephalic Intubation, Canine Jugular Access and Canine Cephalic Venipuncture models
- High quality materials allow for repeated usage and extensive training sessions
- 145x175x165mm Canine Brachycephalic model
- 170x370x345mm Canine Jugular Access model
- 70x100x300mm Canine Cephalic Venipuncture model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Vet Tech Kit	MD-E-C7000S	\$3,307.50

Canine Emergency Neck Access Kit



General Specifications

- Kit includes Canine Cricothyrotomy and Tracheostomy, Canine Jugular Access
- High quality materials allow for repeated usage and extensive training sessions
- 170x370x345mm Canine Cricothyrotomy and Tracheostomy and Canine Jugular Access

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Emergency Neck Access Kit	MD-E-C7001S	\$1,984.00

Canine Airway Management Kit



General Specifications

- Kit includes Canine Brachycephalic Intubation, Canine Cricothyrotomy and Tracheostomy, and Canine Bronchoscopy Trainer
- High quality materials allow for repeated usage and extensive training sessions
- 145x175x165mm Canine Brachycephalic model
- 170x370x345mm Canine Cricothyrotomy and Tracheostomy
- 100x110x470mm Medium bronchoscopy trainer

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Canine Airway Management Kit	MD-E-C7011S	\$2,887.50

Exotic Intubation Kit



General Specifications

- Kit includes Cockatoo Intubation, Macaw Intubation and Rabbit Intubation
- High quality materials allow for repeated usage and extensive training sessions
- 260x130x150mm Cockatoo Intubation model
- 265x155x150mm Macaw Intubation model
- 160x130x75mm Rabbit Intubation model

<u>Component</u>	<u>SKU</u>	<u>Price</u>
Exotic Intubation Kit	MD-E-E7021S	\$1,050.00

Models available for PreOrder*

<u>Products</u>
Canine TPLO
Canine Epidural Standard
Canine Lateral Saphenous Venipuncture
Canine Bandage of the Forelimb Economy
Canine Intraosseous Catheterization Economy
Canine Female Urinary Catheterization Standard
Canine Tooth Extraction Economy
Canine Dental Cleaning

Models in Development**

<u>Products</u>
Canine Intubation Standard w/ NG
Canine Epidural Standard
Canine Saphenous Venipuncture
Feline Saphenous Venipuncture
Blind Cystocentesis
Ultrasound Guided Cystocentesis
Canine Bandage of the Forelimb
Canine Stomach Standard
Canine Intraosseous Catheterization
Canine Joint Injection Thoracic Limb
Canine TPLO
Canine Spay Economy
Canine Tooth Extraction Standard
Canine Dental Cleaning

*PreOrder Models are available for a 10% discount

**Please note this list is ever evolving and therefore non-binding, To suggest a new model or join a focus group, please email info@med-dimensions.com

Index: Vet Bones

<u>Product</u>	<u>SKU</u>	<u>Cost</u>
Canine Humerus Left	MD-E-C0014L	\$27.30
Canine Humerus Right	MD-E-C0014R	\$27.30
Canine Radius Left	MD-E-C0024L	\$27.30
Canine Radius Right	MD-E-C0024R	\$27.30
Canine Ulna Left	MD-E-C0034L	\$25.20
Canine Ulna Right	MD-E-C0034R	\$25.20
Canine Radius/Ulna Left	MD-E-C0044L	\$27.30
Canine Radius/Ulna Right	MD-E-C0044R	\$27.30
Canine Scapula Left	MD-E-C0054L	\$27.30
Canine Scapula Right	MD-E-C0054R	\$27.30
Canine Femur Left	MD-E-C0015L	\$27.30
Canine Femur Right	MD-E-C0015R	\$27.30
Canine Tibia Left	MD-E-C0025L	\$27.30
Canine Tibia Right	MD-E-C0025R	\$27.30
Canine Fibula Left	MD-E-C0035L	\$23.10
Canine Fibula Right	MD-E-C0035R	\$23.10
Canine Tibia/Fibula Left	MD-E-C0045L	\$27.30
Canine Tibia/Fibula Right	MD-E-C0045R	\$27.30
Canine Skull Clear	MD-E-C0011Clear	\$225.75

*Please note prices are subject to change based on material availability and cost. To receive a guaranteed quote, email sales@med-dimensions.com

Index

<u>Product</u>	<u>SKU</u>	<u>Cost</u>
Canine Stifle Visual Aid	MD-E-C0005E	\$133.35
Canine Stifle Visual Aid with Color	MD-E-C0005EColor	\$147.00
3 Slit Suture Pad Economy (Bundle of 5)	MD-E-C1000E	\$96.60
5 Slit Suture Pad Economy(Bundle of 5)	MD-E-C1001E	\$107.10
3 Slit Suture Pad Standard (Bundle of 5)	MD-E-C1000S	\$133.35
5 Slit Suture Pad Standard (Bundle of 5)	MD-E-C1001S	\$144.90
Canine Ear	MD-E-C0011S	\$472.50
Replacement Canine Ear	MD-E-C0011S-1	\$393.75
Feline Intubation Economy	MD-E-F1021E	\$157.50
Feline Intubation Standard	MD-E-F1021S	\$630.00
Replacement Feline Intubation Jaw/Tongue	MD-E-F1021S-1	\$126.00
Replacement Feline Intubation Teeth	MD-E-F1021S-2	\$26.25
Feline Intubation Standard w/ NG	MD-E-F1031S	\$735.00
Replacement Feline Intubation Nose Puck	MD-E-F1031S-1	\$68.25
Feline Intubation Standard Bisected	MD-E-F1041S	\$630.00
Canine Intubation Standard	MD-E-C1021S	\$945.00
Replacement Canine Intubation Jaw/Tongue	MD-E-C1021S-1	\$189.00
Canine Brachycephalic Intubation	MD-E-C1031S	\$997.50
Replacement Brachycephalic Jaw/Tongue	MD-E-C1031S-1	\$194.25
Replacement Brachycephalic Soft Palate	MD-E-C1031S-2	\$42.00
Cockatoo Intubation Standard	MD-E-A1021S	\$367.50

Index

<u>Product</u>	<u>SKU</u>	<u>Cost</u>
Macaw Intubation Standard	MD-E-A1031S	\$367.50
Rabbit Intubation Trainer	MD-E-R1021S	\$393.75
Canine Cephalic Venipuncture	MD-E-C1044S	\$1,155.00
Replacement Cephalic Venipuncture Skin	MD-E-C1044S-1	\$105.00
Replacement Cephalic Vein Beginner 5 Count	MD-E-C1044S-2	\$131.25
Replacement Cephalic Vein Standard 5 Count	MD-E-C1044S-3	\$131.25
Canine Jugular Access	MD-E-C1051S	\$1,260.00
Canine Jugular Central Line Cradle	MD-E-C1051S-1	\$73.50
Replacement Canine Jugular Skin	MD-E-C1051S-2	\$126.00
Replacement Jugular Vessels- 3 Pair	MD-E-C1051S-3	\$147.00
Canine Jugular Access Conversion Kit for Cric/Trache	MD-E-C1051S-4	\$840.00
Canine Cricothyrotomy & Tracheostomy	MD-E-C2031S	\$1,260.00
Replacement Canine Cric & Trache Membrane Tape (12 Count)	MD-E-C2031S-1	\$18.90
Canine Cricothyrotomy & Tracheostomy Conversion Kit for Jugular Access	MD-E-C2031S-3	\$840.00
Canine Bronchoscopy Medium	MD-E-2002SM	\$735.00
Canine Bronchoscopy Large	MD-E-2002SL	\$840.00
Canine Pericardiocentesis Econ	MD-E-C2012E	\$630.00
Canine Pericardiocentesis Standard	MD-E-C2012S	\$787.50
Canine Thoracocentesis	MD-E-C2023	\$2,835.00

Index

<u>Product</u>	<u>SKU</u>	<u>Cost</u>
Canine Joint Injection Pelvic Limb	MD-E-C3005	\$2100.00
Replacement Canine Joint Injection Skin	MD-E-C3005-1	\$157.50
Replacement Canine Joint Injection Pucks	MD-E-C3005-2	\$110.25
Canine Arthroscopy Stifle	MD-E-C3015	\$1,260.00
Replacement Arthroscopy Skin	MD-E-C3015-1	\$157.50
Replacement Arthroscopy Puck	MD-E-C3015-2	\$110.25
Canine Vet Tech Kit	MD-E-C7000S	\$3,307.50
Canine Emergency Neck Access Kit	MD-E-C-7001S	\$1,984.00
Canine Airway Management Kit	MD-E-C7011S	\$2,887.50
Exotic Intubation Kit	MD-E-E7021S	\$1,050.00

Best Practices

To keep your models in best working order we recommend the following

1. Store your model in a cool, dark, and low dust location
2. If you are using any liquid with your model, make sure to wipe it down, leave all connectors open and allow to dry completely before storing
3. Re-apply Silicon powder as needed to reduce tack and dust accumulation
4. Work on a flat surface to prevent the model from slipping

